

## General Device Features Listing and Specifications

<b>External Dimensions</b>	54 cm width x 50 cm depth x 40 cm height
<b>Printing Technologies</b>	Syringe and/or Filament
<b>Number of Tool/Material Bays</b>	1,2, 3, or 4
<b>Primary Fields of Use</b>	Material Science, Tissue Engineering, 3D Cell Culture, Food Science, General Prototyping with software packages for each, as appropriate
<b>Exemplary Materials</b>	plastics, ceramics, silicones, food pastes, hydrogels, organic material, and novel research materials
<b>Usability</b>	Point and click software offers access to fundamental printing process at 'research level' without need for advanced coding or engineering knowledge. Unit is pre-assembled with no training required to operate. Uses standard wall outlet. Disposable and sterilized cartridges available, sterilizable work surface.
<b>Positioning Accuracy</b>	10 µm
<b>Maximum travel speed</b>	130 mm/s
<b>Typical travel speed</b>	80 mm/s
<b>Build Dimensions (x/y/z)**</b>	127mm, 200mm, 65mm
<b>Reservoir volume</b>	Nordson EFD or Becton Dickenson Syringes 3-55 mL
<b>Minimum tip diameter</b>	0.004" / 0.1 mm / 32 gauge (or any lure lock tip)
<b>Maximum tip diameter</b>	0.06" / 1.54 mm / 14 gauge (or any lure lock tip)
<b>File Types</b>	STL and XDFL

\*Specifications may vary based on your unit's specific configuration.



## Tools and Accessories Listing and Specifications

<b>Syringe Tools</b>	Use 3-55mL Nordson EFD or 10ml Becton Dickinson syringes, needles, or taper tips to extrude material using pressure drive. 1,2, 3, or 4 syringe capacity tools available.
<b>Plastic Filament Tool</b>	traditional plastic 3d printing tool (1.75mm or 3mm filament options); Bowden Drive; 80 - 260C range
<b>3D Cell Culture/Well Plating</b>	use a specialized tool head and point-and-click software to easily run, design, or share complex cell culture protocols for automatic 3D cell culture in standard well plates without CAD or complex software/procedures.
<b>Heated build tray</b>	External control (150 C max); useful in plastic printing, cell temperature regulation, and collagen crosslinking.
<b>Syringe Head Heater</b>	Regulate the syringe temperature (80 C max)
<b>Cooled build tray</b>	Peltier junction cooling system (-3.6 C min, 18.4C max)
<b>Syringe Head Cooler</b>	Regulate the syringe temperature (-3.6 C min, 18.4C max)
<b>(UV) LED light tools</b>	LED light source of selected source (e.g. 365 or 385nm), mountable in several configurations. Useful for hydrogel cross linking, among other uses. Illumination coordinated with printing process automatically or manually, depending on tool selected.
<b>USB microscope tool</b>	2 MP camera with 10-40x zoom; can be positioned to view work surface or tool head during printing. (May use multiple per printer.)



### CUSTOM SOLUTIONS

Custom tools are available upon request for users needing specialized functionality. Help others cite your work by requesting that your custom tool is added to our standard accessories list, allowing others to easily build from your research.